

S-100A7L2 Polyclonal Antibody

Cat No: HR1AP8297

For research use only

Overview

Product Name	S-100A7L2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	11kDa
GeneID?Human?	S100A7L2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	S100A7L2; S100A7B; Protein S100-A7-like 2; S100 calcium-binding protein A7-like 2
Background	similarity:Belongs to the S-100 family.,similarity:Contains 2 EF-hand domains.,