

# Mena Polyclonal Antibody

Cat No: HR1AP7865

For research use only

## Overview

Product Name	Mena Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	67kDa
GeneID?Human?	ENAH
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ENAH; MENA; Protein enabled homolog
Background	<p>enabled homolog (Drosophila)(ENAH) Homo sapiens This gene encodes a member of the enabled/ vasodilator-stimulated phosphoprotein. Members of this gene family are involved in actin-based motility. This protein is involved in regulating the assembly of actin filaments and modulates cell adhesion and motility. Alternate splice variants of this gene have been correlated with tumor invasiveness in certain tissues and these variants may serve as prognostic markers. A pseudogene of this gene is found on chromosome 3. [provided by RefSeq, Sep 2016],</p>