

Tensin-1 (phospho Tyr1326) Polyclonal Antibody

Cat No: HR1AP2642

For research use only

Overview

Product Name	Tensin-1 (phospho Tyr1326) Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	185kDa
GeneID?Human?	TNS1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TNS1; TNS; Tensin-1
Background	tensin 1(TNS1) Homo sapiens The protein encoded by this gene localizes to focal adhesions, regions of the plasma membrane where the cell attaches to the extracellular matrix. This protein crosslinks actin filaments and contains a Src homology 2 (SH2) domain, which is often found in molecules involved in signal transduction. This protein is a substrate of calpain II. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Apr 2015],