

## GCP3 Polyclonal Antibody

Cat No: HR1AP3521

For research use only

### Overview

|                       |  |
|-----------------------|--|
| Product Name          | GCP3 Polyclonal Antibody   |
| Source                | Rabbit   |
| Applications          | WB,IHC-p,ELISA   |
| Species Reactivity    | Human,Mouse  |
| Recommended Dilutions |  |
| Immunogen             |  |
| Species               | Rabbit   |
| Storage               | -20°C/1 year   |
| Isotype               |  |
| Clonality             |  |
| Concentration         | 1 mg/ml  |
| Observed band         | 105kDa   |
| GeneID?Human?         | TUBGCP3  |
| Human Swiss-Prot No.  |  |
| Cellular localization |  |
| Alternative Names     | TUBGCP3; GCP3; Gamma-tubulin complex component 3; GCP-3; hGCP3; Gamma-ring complex protein 104 kDa; h104p; hGrip104; Spindle pole body protein Spc98 homolog; hSpc98   |
| Background            | alternative products:Experimental confirmation may be lacking for some isoforms,function:Gamma-tubulin complex is necessary for microtubule nucleation at the centrosome.,similarity:Belongs to the GCP family.,subunit:Gamma-tubulin complex is composed of gamma-tubulin, GCP2, GCP3, GCP4, GCP5 and GCP6.,tissue specificity:Ubiquitously expressed., |