

Tensin-2 (phospho Tyr483) Polyclonal Antibody

Cat No: HR1AP6445

For research use only

Overview

Product Name	Tensin-2 (phospho Tyr483) Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	152kDa
GeneID?Human?	TENC1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TENC1; KIAA1075; TNS2; Tensin-like C1 domain-containing phosphatase; C1 domain-containing phosphatase and tensin homolog; C1-TEN; Tensin-2
Background	tensin 2(TNS2) Homo sapiens The protein encoded by this gene belongs to the tensin family. Tensin is a focal adhesion molecule that binds to actin filaments and participates in signaling pathways. This protein plays a role in regulating cell migration. Alternative splicing occurs at this locus and three transcript variants encoding three distinct isoforms have been identified. [provided by RefSeq, Jul 2008],