

Synphilin-1 Polyclonal Antibody

Cat No: HR1AP4647

For research use only

Overview

Product Name	Synphilin-1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	100kDa
GeneID?Human?	SNCAIP
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SNCAIP; Synphilin-1; Sph1; Alpha-synuclein-interacting protein
Background	<p>synuclein alpha interacting protein(SNCAIP) Homo sapiens This gene encodes a protein containing several protein-protein interaction domains, including ankyrin-like repeats, a coiled-coil domain, and an ATP/GTP-binding motif. The encoded protein interacts with alpha-synuclein in neuronal tissue and may play a role in the formation of cytoplasmic inclusions and neurodegeneration. A mutation in this gene has been associated with Parkinson's disease. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2015],</p>