

Involucrin Polyclonal Antibody

Cat No: HR1AP3745

For research use only

Overview

Product Name	Involucrin Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	68kDa
GeneID?Human?	IVL
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	IVL; Involucrin
Background	involucrin(IVL) Homo sapiens Involucrin, a component of the keratinocyte crosslinked envelope, is found in the cytoplasm and crosslinked to membrane proteins by transglutaminase. This gene is mapped to 1q21, among calpactin I light chain, trichohyalin, profilaggrin, loricrin, and calyculin. [provided by RefSeq, Jul 2008],