

## ? II Tubulin Polyclonal Antibody

Cat No: HR1AP2288

For research use only

### Overview

Product Name	? II Tubulin Polyclonal Antibody
Source	Rabbit
Applications	WB
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	
Observed band	50kDa
GeneID?Human?	TUBB2A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Tubulin beta-2A chain
Background	tubulin beta 2A class IIa(TUBB2A) Homo sapiens Microtubules, key participants in processes such as mitosis and intracellular transport, are composed of heterodimers of alpha- and beta-tubulins. The protein encoded by this gene is a beta-tubulin. Defects in this gene are associated with complex cortical dysplasia with other brain malformations-5. Two transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, Jul 2015],