

Myosin VA Polyclonal Antibody

Cat No: HR1AP7430

For research use only

Overview

Product Name	Myosin VA Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	220kDa
GeneID?Human?	MYO5A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	MYO5A; MYH12; Unconventional myosin-Va; Dilute myosin heavy chain; non-muscle; Myosin heavy chain 12; Myosin-12; Myoxin
Background	<p>myosin VA(MYO5A) Homo sapiens This gene is one of three myosin V heavy-chain genes, belonging to the myosin gene superfamily. Myosin V is a class of actin-based motor proteins involved in cytoplasmic vesicle transport and anchorage, spindle-pole alignment and mRNA translocation. The protein encoded by this gene is abundant in melanocytes and nerve cells. Mutations in this gene cause Griscelli syndrome type-1 (GS1), Griscelli syndrome type-3 (GS3) and neuroectodermal melanolyosomal disease, or Elejalde disease. Multiple alternatively spliced transcript variants encoding different isoforms have been reported, but the full-length nature of some variants has not been determined. [provided by RefSeq, Dec 2008],</p>