

GCP5 Polyclonal Antibody

Cat No: HR1AP3523

For research use only

Overview

Product Name	GCP5 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	118kDa
GeneID?Human?	TUBGCP5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TUBGCP5; GCP5; KIAA1899; Gamma-tubulin complex component 5; GCP-5
Background	function:Gamma-tubulin complex is necessary for microtubule nucleation at the centrosome.,similarity:Belongs to the GCP family.,subunit:Gamma-tubulin complex is composed of gamma-tubulin, GCP2, GCP3, GCP4, GCP5 and GCP6.,tissue specificity:Widely expressed, with highest levels in heart and skeletal muscle and moderate levels in brain.,