

## GCP6 Polyclonal Antibody

Cat No: HR1AP3524

For research use only

### Overview

Product Name	GCP6 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	200kDa
GeneID?Human?	TUBGCP6
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TUBGCP6; GCP6; KIAA1669; Gamma-tubulin complex component 6; GCP-6
Background	tubulin gamma complex associated protein 6(TUBGCP6) Homo sapiens The protein encoded by this gene is part of a large multisubunit complex required for microtubule nucleation at the centrosome. [provided by RefSeq, Jul 2008],