

# MYH14 Polyclonal Antibody

Cat No: HR1AP3986

For research use only

## Overview

Product Name	MYH14 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	228kDa
GeneID?Human?	MYH14
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	MYH14; KIAA2034; FP17425; Myosin-14; Myosin heavy chain 14; Myosin heavy chain; non-muscle IIc; Non-muscle myosin heavy chain IIc; NMHC II-C
Background	myosin heavy chain 14(MYH14) Homo sapiens This gene encodes a member of the myosin superfamily. The protein represents a conventional non-muscle myosin; it should not be confused with the unconventional myosin-14 (MYO14). Myosins are actin-dependent motor proteins with diverse functions including regulation of cytokinesis, cell motility, and cell polarity. Mutations in this gene result in one form of autosomal dominant hearing impairment. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2011],