

?-actin Polyclonal Antibody

Cat No: HR1AP2261

For research use only

Overview

Product Name	?-actin Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p
Species Reactivity	Human,Mouse,Rat,Rb,Ch,Mk,sheep,X
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	43kDa
GeneID?Human?	ACTB
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ACTB; Actin, cytoplasmic 1; Beta-actin
Background	actin beta(ACTB) Homo sapiens This gene encodes one of six different actin proteins. Actins are highly conserved proteins that are involved in cell motility, structure, and integrity. This actin is a major constituent of the contractile apparatus and one of the two nonmuscle cytoskeletal actins. [provided by RefSeq, Jul 2008],