

HSP60 Polyclonal Antibody

Cat No: HR1AP7613

For research use only

Overview

Product Name	HSP60 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	55kDa
GeneID?Human?	HSPA14
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	HSPA14; HSP60; HSP70L1; Heat shock 70 kDa protein 14; HSP70-like protein 1; Heat shock protein HSP60
Background	similarity:Belongs to the heat shock protein 70 family.,