

Arrdc1 Polyclonal Antibody

Cat No: HR1AP2819

For research use only

Overview

Product Name	Arrdc1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	48kDa
GeneID?Human?	ARRDC1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ARRDC1; Arrestin domain-containing protein 1
Background	similarity:Belongs to the arrestin family.,