

DACA-1 Polyclonal Antibody

Cat No: HR1AP3251

For research use only

Overview

Product Name	DACA-1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	61kDa
GeneID?Human?	YTHDF1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	YTHDF1; C20orf21; YTH domain family protein 1; Dermatomyositis associated with cancer putative autoantigen 1; DACA-1
Background	similarity:Contains 1 YTH domain.,