

GAS41 Polyclonal Antibody

Cat No: HR1AP8739

For research use only

Overview

Product Name	GAS41 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	40kDa
GeneID?Human?	YEATS4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	YEATS4; GAS41; YEATS domain-containing protein 4; Glioma-amplified sequence 41; Gas41; NuMA-binding protein 1; NuBI-1; NuBI1
Background	YEATS domain containing 4(YEATS4) Homo sapiens The protein encoded by this gene is found in the nucleoli. It has high sequence homology to human MLLT1, and yeast and human MLLT3 proteins. Both MLLT1 and MLLT3 proteins belong to a class of transcription factors, indicating that the encoded protein might also represent a transcription factor. This protein is thought to be required for RNA transcription. This gene has been shown to be amplified in tumors. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2014].