

AIM2 Polyclonal Antibody

Cat No: HR1AP2728

For research use only

Overview

Product Name	AIM2 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	39kDa
GeneID?Human?	AIM2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	AIM2; Interferon-inducible protein AIM2; Absent in melanoma 2
Background	absent in melanoma 2(AIM2) Homo sapiens AIM2 is a member of the IFI20X /IFI16 family. It plays a putative role in tumorigenic reversion and may control cell proliferation. Interferon-gamma induces expression of AIM2. [provided by RefSeq, Jul 2008],