

IRS-1 (phospho Ser307) Polyclonal Antibody

Cat No: HR1AP2451

For research use only

Overview

| | |
|-----------------------|--|
| Product Name | IRS-1 (phospho Ser307) Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,IHC-p,ELISA |
| Species Reactivity | Human,Mouse,Rat,Monkey |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 170kDa |
| GeneID?Human? | IRS1 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | IRS1; Insulin receptor substrate 1; IRS-1 |
| Background | insulin receptor substrate 1(IRS1) Homo sapiens This gene encodes a protein which is phosphorylated by insulin receptor tyrosine kinase. Mutations in this gene are associated with type II diabetes and susceptibility to insulin resistance. [provided by RefSeq, Nov 2009], |