

Epsin 2 Polyclonal Antibody

Cat No: HR1AP3396

For research use only

Overview

Product Name	Epsin 2 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	68kDa
GeneID?Human?	EPN2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	EPN2; KIAA1065; Epsin-2; EPS-15-interacting protein 2
Background	epsin 2(EPN2) Homo sapiens This gene encodes a protein which interacts with clathrin and adaptor-related protein complex 2, alpha 1 subunit. The protein is found in a brain-derived clathrin-coated vesicle fraction and localizes to the peri-Golgi region and the cell periphery. The protein is thought to be involved in clathrin-mediated endocytosis. Alternate splicing of this gene results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008],