

Melan-A Polyclonal Antibody

Cat No: HR1AP6422

For research use only

Overview

Product Name	Melan-A Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	15kDa
GeneID?Human?	MLANA
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	MLANA; MART1; Melanoma antigen recognized by T-cells 1; MART-1; Antigen LB39-AA; Antigen SK29-AA; Protein Melan-A
Background	tissue specificity:Expression is restricted to melanoma and melanocyte cell lines and retina.,