

# TRIM72 Polyclonal Antibody

Cat No: HR1AP9291

For research use only

## Overview

Product Name	TRIM72 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	
Observed band	53kDa
GeneID?Human?	TRIM72
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Tripartite motif-containing protein 72; Mitsugumin-53; Mg53
Background	similarity:Belongs to the TRIM/RBCC family.,similarity:Contains 1 B box-type zinc finger.,similarity:Contains 1 B30.2/SPRY domain.,similarity:Contains 1 RING-type zinc finger.,