

NUDC (phospho Ser326) Polyclonal Antibody

Cat No: HR1AP5653

For research use only

Overview

Product Name	NUDC (phospho Ser326) Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	45kDa
GeneID?Human?	NUDC
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	NUDC; Nuclear migration protein nudC; Nuclear distribution protein C homolog
Background	nuclear distribution C, dynein complex regulator(NUDC) Homo sapiens This gene encodes a nuclear distribution protein that plays an essential role in mitosis and cytokinesis. The encoded protein is involved in spindle formation during mitosis and in microtubule organization during cytokinesis. Pseudogenes of this gene are found on chromosome 2. [provided by RefSeq, Feb 2012],