

# ANKRD30A Polyclonal Antibody

Cat No: HR1AP2764

For research use only

## Overview

Product Name	ANKRD30A Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	153kDa
GeneID?Human?	ANKRD30A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ANKRD30A; Ankyrin repeat domain-containing protein 30A; Serologically defined breast cancer antigen NY-BR-1
Background	similarity:Contains 6 ANK repeats.,tissue specificity:Mainly expressed in breast and testis. A very faint signal is detected in placenta. Also expressed in many breast cancer cells.,