

Glypican-5 Polyclonal Antibody

Cat No: HR1AP6361

For research use only

Overview

Product Name	Glypican-5 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	65kDa
GeneID?Human?	GPC5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GPC5; Glypican-5
Background	glypican 5(GPC5) Homo sapiens Cell surface heparan sulfate proteoglycans are composed of a membrane-associated protein core substituted with a variable number of heparan sulfate chains. Members of the glypican-related integral membrane proteoglycan family (GRIPS) contain a core protein anchored to the cytoplasmic membrane via a glycosyl phosphatidylinositol linkage. These proteins may play a role in the control of cell division and growth regulation. [provided by RefSeq, Jul 2008].