

## ZnT-2 Polyclonal Antibody

Cat No: HR1AP8662

For research use only

### Overview

Product Name	ZnT-2 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	SLC30A2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SLC30A2; ZNT2; Zinc transporter 2; ZnT-2; Solute carrier family 30 member 2
Background	solute carrier family 30 member 2(SLC30A2) Homo sapiens The protein encoded by this gene is a zinc transporter that acts as a homodimer. The encoded protein plays a role in secreting zinc into breast milk. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2015],