

# ZnT-1 Polyclonal Antibody

Cat No: HR1AP8661

For research use only

## Overview

Product Name	ZnT-1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	55kDa
GeneID?Human?	SLC30A1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SLC30A1; ZNT1; Zinc transporter 1; ZnT-1; Solute carrier family 30 member 1
Background	function:May be involved in zinc transport out of the cell.,similarity:Belongs to the cation diffusion facilitator (CDF) transporter (TC 2.A.4) family. SLC30A subfamily.,subunit:Multimer.,