

TTC23 Polyclonal Antibody

Cat No: HR1AP8316

For research use only

Overview

Product Name	TTC23 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	50kDa
GeneID?Human?	TTC23
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TTC23; HCC8; Tetratricopeptide repeat protein 23; TPR repeat protein 23; Cervical cancer proto-oncogene 8 protein; HCC-8
Background	similarity:Contains 4 TPR repeats.,