

VPS54 Polyclonal Antibody

Cat No: HR1AP7640

For research use only

Overview

| | |
|-----------------------|---|
| Product Name | VPS54 Polyclonal Antibody |
| Source | Rabbit |
| Applications | IHC-p,ELISA |
| Species Reactivity | Human |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | kDa |
| GeneID?Human? | VPS54 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | VPS54; HCC8; Vacuolar protein sorting-associated protein 54; Hepatocellular carcinoma protein 8; Tumor antigen HOM-HCC-8; Tumor antigen SLP-8p |
| Background | VPS54, GARP complex subunit(VPS54) Homo sapiens This gene encodes for a protein that in yeast forms part of a trimeric vacuolar-protein-sorting complex that is required for retrograde transport of proteins from prevacuoles to the late Golgi compartment. As in yeast, mammalian Vps54 proteins contain a coiled-coil region and dileucine motifs. Alternative splicing results in multiple transcript variants encoding different isoforms [provided by RefSeq, Jul 2008]. |