

SVOP Polyclonal Antibody

Cat No: HR1AP7841

For research use only

Overview

Product Name	SVOP Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	61kDa
GeneID?Human?	SVOP
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SVOP; Synaptic vesicle 2-related protein; SV2-related protein
Background	similarity:Belongs to the major facilitator superfamily.,