

## SLC10A7 Polyclonal Antibody

Cat No: HR1AP8836

For research use only

### Overview

Product Name	SLC10A7 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	39kDa
GeneID?Human?	SLC10A7
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SLC10A7; C4orf13; P7; PSEC0051; Sodium/bile acid cotransporter 7; Na(+)/bile acid cotransporter 7; Solute carrier family 10 member 7
Background	similarity:Belongs to the sodium:bile acid symporter family.,