

SLC10A4 Polyclonal Antibody

Cat No: HR1AP6229

For research use only

Overview

Product Name	SLC10A4 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	SLC10A4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SLC10A4; Sodium/bile acid cotransporter 4; Na(+)/bile acid cotransporter 4; Solute carrier family 10 member 4
Background	similarity:Belongs to the sodium:bile acid symporter family.,