

SH3RF3 Polyclonal Antibody

Cat No: HR1AP6944

For research use only

Overview

Product Name	SH3RF3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	92kDa
GeneID?Human?	SH3RF3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SH3RF3; POSH2; SH3MD4; SH3 domain-containing RING finger protein 3; Plenty of SH3s 2; SH3 multiple domains protein 4
Background	domain:The third SH3 domain mediates interaction with PAK2.,similarity:Belongs to the SH3RF family.,similarity:Contains 1 RING-type zinc finger.,similarity:Contains 4 SH3 domains.,subunit:Interacts with PAK2.,