

Serpina6 Polyclonal Antibody

Cat No: HR1AP8957

For research use only

Overview

Product Name	Serpina6 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	SERPINA6
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SERPINA6; CBG; Corticosteroid-binding globulin; CBG; Serpin A6; Transcortin
Background	serpin family A member 6(SERPINA6) Homo sapiens This gene encodes an alpha-globulin protein with corticosteroid-binding properties. This is the major transport protein for glucorticoids and progestins in the blood of most vertebrates. The gene localizes to a chromosomal region containing several closely related serine protease inhibitors which may have evolved by duplication events. [provided by RefSeq, Jul 2008],