

Serpin A11 Polyclonal Antibody

Cat No: HR1AP6544

For research use only

Overview

Product Name	Serpin A11 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	47kDa
GeneID?Human?	SERPINA11
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SERPINA11; Serpin A11
Background	similarity:Belongs to the serpin family.,