

Rsu-1 Polyclonal Antibody

Cat No: HR1AP8221

For research use only

Overview

Product Name	Rsu-1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	36kDa
GeneID?Human?	RSU1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RSU1; RSP1; Ras suppressor protein 1; RSP-1; Rsu-1
Background	<p>Ras suppressor protein 1(RSU1) Homo sapiens This gene encodes a protein that is involved in the Ras signal transduction pathway, growth inhibition, and nerve-growth factor induced differentiation processes, as determined in mouse and human cell line studies. In mouse, the encoded protein was initially isolated based on its ability to inhibit v-Ras transformation. Multiple alternatively spliced transcript variants for this gene have been reported; one of these variants was found only in glioma tumors. [provided by RefSeq, Jul 2008].</p>