

# Ribosomal Protein L7L1 Polyclonal Antibody

Cat No: HR1AP6724

For research use only

## Overview

Product Name	Ribosomal Protein L7L1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	29kDa
GeneID?Human?	RPL7L1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RPL7L1; 60S ribosomal protein L7-like 1
Background	similarity:Belongs to the ribosomal protein L30P family.,