

# Ribosomal Protein L3L Polyclonal Antibody

Cat No: HR1AP8174

For research use only

## Overview

Product Name	Ribosomal Protein L3L Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	RPL3L
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RPL3L; 60S ribosomal protein L3-like
Background	<p>ribosomal protein L3 like(RPL3L) Homo sapiens This gene encodes a protein that shares sequence similarity with ribosomal protein L3. The protein belongs to the L3P family of ribosomal proteins. Unlike the ubiquitous expression of ribosomal protein genes, this gene has a tissue-specific pattern of expression, with the highest levels of expression in skeletal muscle and heart. It is not currently known whether the encoded protein is a functional ribosomal protein or whether it has evolved a function that is independent of the ribosome. [provided by RefSeq, Jul 2008].</p>