

Ribosomal Protein L39L Polyclonal Antibody

Cat No: HR1AP5774

For research use only

Overview

Product Name	Ribosomal Protein L39L Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	RPL39L
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RPL39L; RPL39L1; 60S ribosomal protein L39-like; 60S ribosomal protein L39-2
Background	ribosomal protein L39 like(RPL39L) Homo sapiens This gene encodes a protein sharing high sequence similarity with ribosomal protein L39. Although the name of this gene has been referred to as 'ribosomal protein L39' in the public databases, its official name is 'ribosomal protein L39-like'. It is not currently known whether the encoded protein is a functional ribosomal protein or whether it has evolved a function that is independent of the ribosome. [provided by RefSeq, Jul 2008],