

Ribosomal Protein L37A Polyclonal Antibody

Cat No: HR1AP8184

For research use only

Overview

Product Name	Ribosomal Protein L37A Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	10kDa
GeneID?Human?	RPL37A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RPL37A; 60S ribosomal protein L37a
Background	<p>ribosomal protein L37a(RPL37A) Homo sapiens Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 60S subunit. The protein belongs to the L37AE family of ribosomal proteins. It is located in the cytoplasm. The protein contains a C4-type zinc finger-like domain. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome. [provided by RefSeq, Jul 2008],</p>