

## RB3 Polyclonal Antibody

Cat No: HR1AP8755

For research use only

### Overview

Product Name	RB3 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	22kDa
GeneID?Human?	STMN4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	STMN4; Stathmin-4; Stathmin-like protein B3; RB3
Background	similarity:Belongs to the stathmin family.,