

RFPL4B Polyclonal Antibody

Cat No: HR1AP7391

For research use only

Overview

Product Name	RFPL4B Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	RFPL4B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RFPL4B; RNF211; Ret finger protein-like 4B; RING finger protein 211
Background	similarity:Contains 1 B30.2/SPRY domain.,similarity:Contains 1 RING-type zinc finger.,