

REPS1 Polyclonal Antibody

Cat No: HR1AP8910

For research use only

Overview

Product Name	REPS1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	80kDa
GeneID?Human?	REPS1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	REPS1; RalBP1-associated Eps domain-containing protein 1; RalBP1-interacting protein 1
Background	<p>RALBP1 associated Eps domain containing 1(REPS1) Homo sapiens This gene encodes a signaling adaptor protein with two EH domains that interacts with proteins that participate in signaling, endocytosis and cytoskeletal changes. The encoded protein has been found in association with intersectin 1 and Src homology 3-domain growth factor receptor-bound 2-like (endophilin) interacting protein 1 when intersectin 1 was isolated from clathrin-coated pits. The encoded protein has also been shown to interact with amphiphysin, a cytoplasmic protein at the surface of synaptic vesicles. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Mar 2014].</p>