

Pxmp2 Polyclonal Antibody

Cat No: HR1AP8062

For research use only

Overview

Product Name	Pxmp2 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	PXMP2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	PXMP2; PMP22; Peroxisomal membrane protein 2; 22 kDa peroxisomal membrane protein
Background	function:Seems to be involved in pore-forming activity and may contribute to the unspecific permeability of the peroxisomal membrane.,similarity:Belongs to the peroxisomal membrane protein PXMP2/4 family.,subunit:Interacts with PEX19 and SIVA1.,