

Prefoldin 1 Polyclonal Antibody

Cat No: HR1AP7685

For research use only

Overview

Product Name	Prefoldin 1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	14kDa
GeneID?Human?	PFDN1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	PFDN1; PFD1; Prefoldin subunit 1
Background	<p>prefoldin subunit 1(PFDN1) Homo sapiens This gene encodes a member of the prefoldin beta subunit family. The encoded protein is one of six subunits of prefoldin, a molecular chaperone complex that binds and stabilizes newly synthesized polypeptides, thereby allowing them to fold correctly. The complex, consisting of two alpha and four beta subunits, forms a double beta barrel assembly with six protruding coiled-coils. [provided by RefSeq, Jul 2008],</p>