

Podocan Polyclonal Antibody

Cat No: HR1AP5839

For research use only

Overview

Product Name	Podocan Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	PODN
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	PODN; SLRR5A; Podocan
Background	podocan(PODN) Homo sapiens The protein encoded by this gene is a member of the small leucine-rich repeat protein family and contains an amino terminal CX3CX7C cysteine-rich cluster followed by a leucine-rich repeat domain. Studies suggest that this protein could function to inhibit smooth muscle cell proliferation and migration following arterial injury. [provided by RefSeq, Jul 2016],