

# Phocein Polyclonal Antibody

Cat No: HR1AP6551

For research use only

## Overview

|                       |   |
|-----------------------|---|
| Product Name          | Phocein Polyclonal Antibody   |
| Source                | Rabbit  |
| Applications          | WB,IHC-p,IF,ELISA   |
| Species Reactivity    | Human,Mouse,Rat   |
| Recommended Dilutions |   |
| Immunogen             |   |
| Species               | Rabbit  |
| Storage               | -20°C/1 year  |
| Isotype               |   |
| Clonality             |   |
| Concentration         | 1 mg/ml   |
| Observed band         | 26kDa   |
| GeneID?Human?         | MOB4  |
| Human Swiss-Prot No.  |   |
| Cellular localization |   |
| Alternative Names     | MOB4; MOB3; MOBKL3; PHOCN; PREI3; CGI-95; MOB-like protein phocein; 2C4D; Class II mMOB1; Mob1 homolog 3; Mob3; Mps one binder kinase activator-like 3; Preimplantation protein 3   |
| Background            | MOB family member 4, phocein(MOB4) Homo sapiens This gene was identified based on its similarity with the mouse counterpart. Studies of the mouse counterpart suggest that the expression of this gene may be regulated during oocyte maturation and preimplantation following zygotic gene activation. Alternatively spliced transcript variants encoding distinct isoforms have been observed. Naturally occurring read-through transcription occurs between this locus and the neighboring locus HSPE1.[provided by RefSeq, Feb 2011], |