

## Pki ? Polyclonal Antibody

Cat No: HR1AP7856

For research use only

### Overview

Product Name	Pki ? Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,IF,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	PKIA
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	PKIA; PRKACN1; cAMP-dependent protein kinase inhibitor alpha; PKI-alpha; cAMP-dependent protein kinase inhibitor; muscle/brain isoform
Background	protein kinase (cAMP-dependent, catalytic) inhibitor alpha(PKIA) Homo sapiens The protein encoded by this gene is a member of the cAMP-dependent protein kinase (PKA) inhibitor family. This protein was demonstrated to interact with and inhibit the activities of both C alpha and C beta catalytic subunits of the PKA. Alternatively spliced transcript variants encoding the same protein have been reported. [provided by RefSeq, Jul 2008],