

NDRG3 Polyclonal Antibody

Cat No: HR1AP8006

For research use only

Overview

Product Name	NDRG3 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	45kDa
GeneID?Human?	NDRG3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	NDRG3; Protein NDRG3; N-myc downstream-regulated gene 3 protein
Background	similarity:Belongs to the NDRG family.,tissue specificity:Ubiquitous. Highly expressed in brain.,