

# NDRG4 Polyclonal Antibody

Cat No: HR1AP8322

For research use only

## Overview

Product Name	NDRG4 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	34kDa
GeneID?Human?	NDRG4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	NDRG4; BDM1; KIAA1180; Protein NDRG4; Brain development-related molecule 1; N-myc downstream-regulated gene 4 protein; Vascular smooth muscle cell-associated protein 8; SMAP-8
Background	NDRG family member 4(NDRG4) Homo sapiens This gene is a member of the N-myc downregulated gene family which belongs to the alpha/beta hydrolase superfamily. The protein encoded by this gene is a cytoplasmic protein that is required for cell cycle progression and survival in primary astrocytes and may be involved in the regulation of mitogenic signalling in vascular smooth muscles cells. Alternative splicing results in multiple transcripts encoding different isoforms.[provided by RefSeq, Jun 2011],